

# Product Catalogue

Medical Compression Stockings and Accessories



Comfort, Health & Style!



# **JOBST<sup>®</sup> UltraSheer**

Medical Compression Stockings



**The exceptional  
combination of  
comfort, style  
and gradient  
compression  
therapy**



# JOBST® UltraSheer

## Medical Compression Stockings

### Sizing Chart:

Size	Ankle 	Calf 	Thigh 	Hip 	Hip (Maternity)
<b>S</b>	18-21 cm	28-38 cm	40-62 cm	71-117 cm	81-132 cm
<b>M</b>	21-25 cm	30-42 cm	46-70 cm	76-127 cm	86-142 cm
<b>L</b>	25-29 cm	32-46 cm	54-78 cm	81-137 cm	91-152 cm
<b>XL</b>	29-33 cm	34-50 cm	60-81 cm	102-166 cm	112-180 cm
<b>LFC</b>	25-30 cm	46-61 cm	–	–	–
<b>XLFC</b>	30-36 cm	46-61 cm	–	–	–

Regular : ≥ 38 cm	Petite : < 38 cm	
Regular : ≥ 68.5 cm	Petite : < 68.5 cm	

For Knee High, measure Ankle and Calf. For Thigh High, measure Ankle, Calf and Thigh.  
For Pantyhose and Maternity, measure Ankle, Thigh and Hip.

### Available Colors



Natural



Honey



Sun  
Bronze



Suntan



Black



Espresso



Anthracite

### Petite Knee High (Closed Toe)



Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7540802	–	–	7540803	7540801	–	–
15-20 mmHg	<b>M</b>	7540902	–	–	7540903	7540901	–	–
15-20 mmHg	<b>L</b>	7541002	–	–	7541003	7541001	–	–
15-20 mmHg	<b>XL</b>	7541502	–	–	7541503	7541501	–	–
20-30 mmHg	<b>S</b>	7543502	–	–	7543503	7543501	–	–
20-30 mmHg	<b>M</b>	7543602	–	–	7543603	7543601	–	–
20-30 mmHg	<b>L</b>	7543702	–	–	7543703	7543701	–	–
20-30 mmHg	<b>XL</b>	7543802	–	–	7543803	7543801	–	–
30-40 mmHg	<b>S</b>	7544902	–	–	7544903	7544901	–	–
30-40 mmHg	<b>M</b>	7545002	–	–	7545003	7545001	–	–
30-40 mmHg	<b>L</b>	7545102	–	–	7545103	7545101	–	–
30-40 mmHg	<b>XL</b>	7545202	–	–	7545203	7545201	–	–

**Fiber Contents:** 15-20 mmHg\*: 75% Nylon, 25% Spandex.  
20-30 mmHg\*: 74% Nylon, 26% Spandex.  
30-40 mmHg\*: 60% Nylon, 40% Spandex.

### Petite Knee High (Open Toe)



Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
20-30 mmHg	<b>S</b>	7543505	–	–	–	7543504	–	–
20-30 mmHg	<b>M</b>	7543605	–	–	–	7543604	–	–
20-30 mmHg	<b>L</b>	7543705	–	–	–	7543704	–	–
20-30 mmHg	<b>XL</b>	7543805	–	–	–	7543804	–	–

**Fiber Contents:** 20-30 mmHg\*: 74% Nylon, 26% Spandex.

# JOBST® UltraSheer

## Medical Compression Stockings



### Knee High (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7516525	7516528	7516531	7516526	7516504	7516529	7516527
15-20 mmHg	<b>M</b>	7516625	7516628	7516631	7516626	7516604	7516629	7516627
15-20 mmHg	<b>L</b>	7516725	7516728	7516731	7516726	7516704	7516729	7516727
15-20 mmHg	<b>XL</b>	7516825	7516828	7516831	7516826	7516804	7516829	7516827
15-20 mmHg	<b>LFC</b>	7541602	–	–	–	7541601	–	–
15-20 mmHg	<b>XLFC</b>	7541702	–	–	–	7541701	–	–
20-30 mmHg	<b>S</b>	7520119	7520123	7520130	7520120	7520104	7520121	7520122
20-30 mmHg	<b>M</b>	7520219	7520223	7520230	7520220	7520204	7520221	7520222
20-30 mmHg	<b>L</b>	7520319	7520323	7520330	7520320	7520304	7520321	7520322
20-30 mmHg	<b>XL</b>	7520419	7520423	7520430	7520420	7520404	7520421	7520422
20-30 mmHg	<b>LFC</b>	7543302	–	–	–	7543301	–	–
20-30 mmHg	<b>XLFC</b>	7543402	–	–	–	7543401	–	–
30-40 mmHg	<b>S</b>	7524121	7524124	7524129	7524122	7524104	7524125	7524123
30-40 mmHg	<b>M</b>	7524221	7524224	7524229	7524222	7524204	7524225	7524223
30-40 mmHg	<b>L</b>	7524321	7524324	7524329	7524322	7524304	7524325	7524323
30-40 mmHg	<b>XL</b>	7524421	7524424	7524429	7524422	7524404	7524425	7524423
30-40 mmHg	<b>LFC</b>	7544702	–	–	–	7544701	–	–
30-40 mmHg	<b>XLFC</b>	7544802	–	–	–	7544801	–	–

**Fiber Contents:** 15-20 mmHg\*: 75% Nylon, 25% Spandex.  
 20-30 mmHg\*: 74% Nylon, 26% Spandex.  
 30-40 mmHg\*: 60% Nylon, 40% Spandex.



### Knee High (Open Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7516530	–	–	–	7516507	–	–
15-20 mmHg	<b>M</b>	7516630	–	–	–	7516607	–	–
15-20 mmHg	<b>L</b>	7516730	–	–	–	7516707	–	–
15-20 mmHg	<b>XL</b>	7516830	–	–	–	7516807	–	–
20-30 mmHg	<b>S</b>	7520129	7520128	–	–	7520127	–	–
20-30 mmHg	<b>M</b>	7520229	7520228	–	–	7520227	–	–
20-30 mmHg	<b>L</b>	7520329	7520328	–	–	7520327	–	–
20-30 mmHg	<b>XL</b>	7520429	7520428	–	–	7520427	–	–
20-30 mmHg	<b>LFC</b>	7543329	7543328	–	–	7543327	–	–
20-30 mmHg	<b>XLFC</b>	7543429	7543428	–	–	7543427	–	–
30-40 mmHg	<b>S</b>	7524128	–	–	–	7524127	–	–
30-40 mmHg	<b>M</b>	7524228	–	–	–	7524227	–	–
30-40 mmHg	<b>L</b>	7524328	–	–	–	7524327	–	–
30-40 mmHg	<b>XL</b>	7524428	–	–	–	7524427	–	–

**Fiber Contents:** 15-20 mmHg\*: 75% Nylon, 25% Spandex.  
 20-30 mmHg\*: 74% Nylon, 26% Spandex.  
 30-40 mmHg\*: 60% Nylon, 40% Spandex.

# JOBST® UltraSheer

## Medical Compression Stockings



### Petite Thigh High (Silicone Lace Band / Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7542902	–	–	7542903	7542901	–	–
15-20 mmHg	<b>M</b>	7543002	–	–	7543003	7543001	–	–
15-20 mmHg	<b>L</b>	7543102	–	–	7543103	7543101	–	–
15-20 mmHg	<b>XL</b>	7543202	–	–	7543203	7543201	–	–
20-30 mmHg	<b>S</b>	7544302	–	–	7544303	7544301	–	–
20-30 mmHg	<b>M</b>	7544402	–	–	7544403	7544401	–	–
20-30 mmHg	<b>L</b>	7544502	–	–	7544503	7544501	–	–
20-30 mmHg	<b>XL</b>	7544602	–	–	7544603	7544601	–	–
30-40 mmHg	<b>S</b>	7545702	–	–	7545703	7545701	–	–
30-40 mmHg	<b>M</b>	7545802	–	–	7545803	7545801	–	–
30-40 mmHg	<b>L</b>	7545902	–	–	7545903	7545901	–	–
30-40 mmHg	<b>XL</b>	7546002	–	–	7546003	7546001	–	–

**Fiber Contents:** 15-20 mmHg\*: 60% Nylon, 25% Spandex, 15% Silicone.  
 20-30 mmHg\*: 59% Nylon, 29% Spandex, 12% Silicone.  
 30-40 mmHg\*: 55% Nylon, 37% Spandex, 8% Silicone.



### Thigh High (Silicone Lace Band / Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7516922	7516925	7516927	7516923	7516904	7516926	7516924
15-20 mmHg	<b>M</b>	7517022	7517025	7517027	7517023	7517004	7517026	7517024
15-20 mmHg	<b>L</b>	7517122	7517125	7517127	7517123	7517104	7517126	7517124
15-20 mmHg	<b>XL</b>	7517222	7517225	7517227	7517223	7517204	7517226	7517224
20-30 mmHg	<b>S</b>	7520523	7520527	7520531	7520524	7520504	7520525	7520526
20-30 mmHg	<b>M</b>	7520623	7520627	7520631	7520624	7520604	7520625	7520626
20-30 mmHg	<b>L</b>	7520723	7520727	7520731	7520724	7520704	7520725	7520726
20-30 mmHg	<b>XL</b>	7520823	7520827	7520831	7520824	7520804	7520825	7520826
30-40 mmHg	<b>S</b>	7524520	7524523	7524525	7524521	7524504	7524524	7524522
30-40 mmHg	<b>M</b>	7524620	7524623	7524625	7524621	7524604	7524624	7524622
30-40 mmHg	<b>L</b>	7524720	7524723	7524725	7524721	7524704	7524724	7524722
30-40 mmHg	<b>XL</b>	7524820	7524823	7524825	7524821	7524804	7524824	7524822

**Fiber Contents:** 15-20 mmHg\*: 60% Nylon, 25% Spandex, 15% Silicone.  
 20-30 mmHg\*: 59% Nylon, 29% Spandex, 12% Silicone.  
 30-40 mmHg\*: 55% Nylon, 37% Spandex, 8% Silicone.

# JOBST® UltraSheer

## Medical Compression Stockings



### Petite Thigh High (Silicone Dot Band / Open Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
20-30 mmHg	<b>S</b>	7543907	–	–	–	–	–	–
20-30 mmHg	<b>M</b>	7544007	–	–	–	–	–	–
20-30 mmHg	<b>L</b>	7544107	–	–	–	–	–	–
20-30 mmHg	<b>XL</b>	7544207	–	–	–	–	–	–

**Fiber Contents:** 20-30 mmHg\*: 64% Nylon, 24% Spandex, 12% Silicone.



### Thigh High (Silicone Dot Band / Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7542502	–	7542506	7542503	7542501	–	–
15-20 mmHg	<b>M</b>	7542602	–	7542606	7542603	7542601	–	–
15-20 mmHg	<b>L</b>	7542702	–	7542706	7542703	7542701	–	–
15-20 mmHg	<b>XL</b>	7542802	–	7542806	7542803	7542801	–	–
20-30 mmHg	<b>S</b>	7543902	–	7543908	7543903	7543901	–	–
20-30 mmHg	<b>M</b>	7544002	–	7544008	7544003	7544001	–	–
20-30 mmHg	<b>L</b>	7544102	–	7544108	7544103	7544101	–	–
20-30 mmHg	<b>XL</b>	7544202	–	7544208	7544203	7544201	–	–
30-40 mmHg	<b>S</b>	7545302	–	7545305	7545303	7545301	–	–
30-40 mmHg	<b>M</b>	7545402	–	7545405	7545403	7545401	–	–
30-40 mmHg	<b>L</b>	7545502	–	7545505	7545503	7545501	–	–
30-40 mmHg	<b>XL</b>	7545602	–	7545605	7545603	7545601	–	–

**Fiber Contents:** 15-20 mmHg\*: 66% Nylon, 19% Spandex, 15% Silicone.

20-30 mmHg\*: 64% Nylon, 24% Spandex, 12% Silicone.

30-40 mmHg\*: 56% Nylon, 34% Spandex, 10% Silicone.



### Thigh High (Silicone Dot Band / Open Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7542505	–	–	–	7542504	–	–
15-20 mmHg	<b>M</b>	7542605	–	–	–	7542604	–	–
15-20 mmHg	<b>L</b>	7542705	–	–	–	7542704	–	–
15-20 mmHg	<b>XL</b>	7542805	–	–	–	7542804	–	–
20-30 mmHg	<b>S</b>	7543906	7543905	–	–	7543904	–	–
20-30 mmHg	<b>M</b>	7544006	7544005	–	–	7544004	–	–
20-30 mmHg	<b>L</b>	7544106	7544105	–	–	7544104	–	–
20-30 mmHg	<b>XL</b>	7544206	7544205	–	–	7544204	–	–
30-40 mmHg	<b>S</b>	7545304	–	–	–	–	–	–
30-40 mmHg	<b>M</b>	7545404	–	–	–	–	–	–
30-40 mmHg	<b>L</b>	7545504	–	–	–	–	–	–
30-40 mmHg	<b>XL</b>	7545604	–	–	–	–	–	–

**Fiber Contents:** 15-20 mmHg\*: 66% Nylon, 19% Spandex, 15% Silicone.

20-30 mmHg\*: 64% Nylon, 24% Spandex, 12% Silicone.

30-40 mmHg\*: 56% Nylon, 34% Spandex, 10% Silicone.

# JOBST® UltraSheer

## Medical Compression Stockings



### Pantyhose (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7517312	7517315	7517317	7517313	7517304	7517316	7517314
15-20 mmHg	<b>M</b>	7517412	7517415	7517417	7517413	7517404	7517416	7517414
15-20 mmHg	<b>L</b>	7517512	7517515	7517517	7517513	7517504	7517516	7517514
15-20 mmHg	<b>XL</b>	7517612	7517615	7517617	7517613	7517604	7517616	7517614
20-30 mmHg	<b>S</b>	7520913	7520917	7520918	7520914	7520904	7520915	7520916
20-30 mmHg	<b>M</b>	7521013	7521017	7521018	7521014	7521004	7521015	7521016
20-30 mmHg	<b>L</b>	7521113	7521117	7521118	7521114	7521104	7521115	7521116
20-30 mmHg	<b>XL</b>	7521213	7521217	7521218	7521214	7521204	7521215	7521216
30-40 mmHg	<b>S</b>	7524913	7524916	7524918	7524914	7524904	7524917	7524915
30-40 mmHg	<b>M</b>	7525013	7525016	7525018	7525014	7525004	7525017	7525015
30-40 mmHg	<b>L</b>	7525113	7525116	7525118	7525114	7525104	7525117	7525115
30-40 mmHg	<b>XL</b>	7525213	7525216	7525218	7525214	7525204	7525217	7525215

**Fiber Contents:** 15-20 mmHg\*: 72% Nylon, 28% Spandex.  
 20-30 mmHg\*: 73% Nylon, 27% Spandex.  
 30-40 mmHg\*: 61% Nylon, 39% Spandex.



### Maternity (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Sun Bronze	Suntan	Black	Espresso	Anthracite
15-20 mmHg	<b>S</b>	7517704	–	–	–	7517701	–	–
15-20 mmHg	<b>M</b>	7517804	–	–	–	7517801	–	–
15-20 mmHg	<b>L</b>	7517904	–	–	–	7517901	–	–
15-20 mmHg	<b>XL</b>	7518004	–	–	–	7518001	–	–
20-30 mmHg	<b>S</b>	7521306	–	–	–	7521303	–	–
20-30 mmHg	<b>M</b>	7521406	–	–	–	7521403	–	–
20-30 mmHg	<b>L</b>	7521506	–	–	–	7521503	–	–
20-30 mmHg	<b>XL</b>	7521606	–	–	–	7521603	–	–

**Fiber Contents:** 15-20 mmHg\*: 72% Nylon, 28% Spandex.  
 20-30 mmHg\*: 76% Nylon, 24% Spandex.

#### Care Instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size.



# JOBST® UltraSheer

with an elegant pattern







The exceptional  
combination of  
comfort, style  
and gradient  
compression  
therapy



# JOBST® UltraSheer

with an elegant pattern

## Sizing chart:

Size	Ankle 	Calf 	Thigh 	Hip 
<b>S</b>	18 - 21 cm	28 - 38 cm	40 - 62 cm	71 - 117 cm
<b>M</b>	21 - 25 cm	30 - 42 cm	46 - 70 cm	76 - 127 cm
<b>L</b>	25 - 29 cm	32 - 46 cm	54 - 78 cm	81 - 137 cm
<b>XL</b>	29 - 33 cm	34 - 50 cm	60 - 81 cm	102 - 166 cm

Regular :  $\geq 38$  cm

Petite :  $< 38$  cm



Regular :  $\geq 68.5$  cm

Petite :  $< 68.5$  cm



For Knee High, measure Ankle and Calf. For Thigh High, measure Ankle, Calf and Thigh.  
For Panty, measure Ankle, Thigh and Hip.

## Available pattern:

Diamond shape 

## Available colors:



Black



Espresso





## Knee High (Closed Toe)

Compression mmHg*	Size	Diamond Black	Diamond Espresso
15 - 20 mmHg	<b>S</b>	7363300	7279600
15 - 20 mmHg	<b>M</b>	7363301	7279601
15 - 20 mmHg	<b>L</b>	7363302	7279602
15 - 20 mmHg	<b>XL</b>	7363303	7279603
20 - 30 mmHg	<b>S</b>	7363400	7279700
20 - 30 mmHg	<b>M</b>	7363401	7279701
20 - 30 mmHg	<b>L</b>	7363402	7279702
20 - 30 mmHg	<b>XL</b>	7363403	7279703

### Fiber contents:

15-20 mmHg\* : 75% Nylon and 25% Spandex.

20-30 mmHg\* : 74% Nylon and 26% Spandex.



## Thigh High (Silicone Dot Band / Closed Toe)

Compression mmHg*	Size	Diamond Black	Diamond Espresso
15 - 20 mmHg	<b>S</b>	7363500	7279604
15 - 20 mmHg	<b>M</b>	7363501	7279605
15 - 20 mmHg	<b>L</b>	7363502	7279606
15 - 20 mmHg	<b>XL</b>	7363503	7279607
20 - 30 mmHg	<b>S</b>	7363600	7279704
20 - 30 mmHg	<b>M</b>	7363601	7279705
20 - 30 mmHg	<b>L</b>	7363602	7279706
20 - 30 mmHg	<b>XL</b>	7363603	7279707

### Fiber contents:

15-20 mmHg\* : 66% Nylon, 19% Spandex, 15% Silicone.

20-30 mmHg\* : 64% Nylon, 24% Spandex, 12% Silicone.

# JOBST® UltraSheer

with an elegant pattern



## Pantyhose (Closed Toe)

Compression mmHg*	Size	Diamond Black	Diamond Espresso
15 - 20 mmHg	<b>S</b>	7363700	7279608
15 - 20 mmHg	<b>M</b>	7363701	7279609
15 - 20 mmHg	<b>L</b>	7363702	7279610
15 - 20 mmHg	<b>XL</b>	7363703	7279611
20 - 30 mmHg	<b>S</b>	7363800	7279708
20 - 30 mmHg	<b>M</b>	7363801	7279709
20 - 30 mmHg	<b>L</b>	7363802	7279710
20 - 30 mmHg	<b>XL</b>	7363803	7279711

### Fiber contents:

15-20 mmHg\*: 72% Nylon and 28% Spandex.

20-30 mmHg\*: 73% Nylon and 27% Spandex.

### Care instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size.

# **JOBST<sup>®</sup> Opaque**

Medical Compression Stockings



**Softness and ease  
of donning  
for exceptional  
comfortable  
leg therapy**



# JOBST® Opaque

## Medical Compression Stockings

### Features & Benefits

- High quality yarns and unique finishes give JOBST® Opaque a perfectly soft touch
- JOBST® Opaque wicks moisture away from the skin effectively and lets it evaporate quickly
- Anatomical designs and silicone blend softeners provide great fit and make JOBST® Opaque easier to don than ever
- JOBST® Opaque's refined opacity discreetly covers blemishes, giving your legs an elegant look
- JOBST® Opaque's durability will provide your patients with long-lasting elegance

### Available Colors



Natural



Honey







Espresso



Black



### Sizing Chart

Size	Ankle 	Calf 	Thigh 	Hip 
<b>S</b>	18 - 21 cm	28 - 38 cm	40 - 62 cm	71 - 117 cm
<b>M</b>	21 - 25 cm	30 - 42 cm	46 - 70 cm	76 - 127 cm
<b>L</b>	25 - 29 cm	32 - 46 cm	54 - 78 cm	81 - 137 cm
<b>XL</b>	29 - 33 cm	34 - 50 cm	60 - 81 cm	102 - 166 cm
<b>LFC*</b>	25 - 30 cm	46 - 61 cm	—	—
<b>XLFC*</b>	30 - 36 cm	46 - 61 cm	—	—

Regular :  $\geq 38$  cm

Petite :  $< 38$  cm



Regular :  $\geq 68.5$  cm

Petite :  $< 68.5$  cm



For Knee High, measure Ankle and Calf.

For Thigh High, measure Ankle, Calf and Thigh.

For Pantyhose, measure Ankle, Thigh and Hip.

# JOBST® Opaque

## Medical Compression Stockings



### Petite Knee High (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518122	–	7518123	–
15-20 mmHg	<b>M</b>	7518222	–	7518223	–
15-20 mmHg	<b>L</b>	7518322	–	7518323	–
15-20 mmHg	<b>XL</b>	7518422	–	7518423	–
20-30 mmHg	<b>S</b>	7521722	–	7521723	–
20-30 mmHg	<b>M</b>	7521822	–	7521823	–
20-30 mmHg	<b>L</b>	7521922	–	7521923	–
20-30 mmHg	<b>XL</b>	7522022	–	7522023	–
30-40 mmHg	<b>S</b>	7522922	–	7522923	–
30-40 mmHg	<b>M</b>	7523022	–	7523023	–
30-40 mmHg	<b>L</b>	7523122	–	7523123	–
30-40 mmHg	<b>XL</b>	7523222	–	7523223	–

#### Fiber Contents:

15-20 mmHg\*: 78% Nylon, 22% Spandex.  
 20-30 mmHg\*: 78% Nylon, 22% Spandex.  
 30-40 mmHg\*: 70% Nylon, 30% Spandex.



### Petite Knee High (Open Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518120	–	–	–
15-20 mmHg	<b>M</b>	7518220	–	–	–
15-20 mmHg	<b>L</b>	7518320	–	–	–
15-20 mmHg	<b>XL</b>	7518420	–	–	–
20-30 mmHg	<b>S</b>	7521720	–	–	–
20-30 mmHg	<b>M</b>	7521820	–	–	–
20-30 mmHg	<b>L</b>	7521920	–	–	–
20-30 mmHg	<b>XL</b>	7522020	–	–	–
30-40 mmHg	<b>S</b>	7522920	–	–	–
30-40 mmHg	<b>M</b>	7523020	–	–	–
30-40 mmHg	<b>L</b>	7523120	–	–	–
30-40 mmHg	<b>XL</b>	7523220	–	–	–

#### Fiber Contents:

15-20 mmHg\*: 78% Nylon, 22% Spandex.  
 20-30 mmHg\*: 78% Nylon, 22% Spandex.  
 30-40 mmHg\*: 70% Nylon, 30% Spandex.

# JOBST® Opaque

## Medical Compression Stockings



### Knee High (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518100	7518109	7518101	7518110
15-20 mmHg	<b>M</b>	7518200	7518209	7518201	7518210
15-20 mmHg	<b>L</b>	7518300	7518309	7518301	7518310
15-20 mmHg	<b>XL</b>	7518400	7518409	7518401	7518410
20-30 mmHg	<b>S</b>	7521700	7521711	7521703	7521713
20-30 mmHg	<b>M</b>	7521800	7521811	7521803	7521813
20-30 mmHg	<b>L</b>	7521900	7521911	7521903	7521913
20-30 mmHg	<b>XL</b>	7522000	7522011	7522003	7522013
20-30 mmHg	<b>LFC</b>	7848922	–	7848923	–
20-30 mmHg	<b>XLFC</b>	7848822	–	7848823	–
30-40 mmHg	<b>S</b>	7522900	7522908	7522902	7522910
30-40 mmHg	<b>M</b>	7523000	7523008	7523002	7523010
30-40 mmHg	<b>L</b>	7523100	7523108	7523102	7523110
30-40 mmHg	<b>XL</b>	7523200	7523208	7523202	7523210
30-40 mmHg	<b>LFC</b>	7848722	–	7848723	–
30-40 mmHg	<b>XLFC</b>	7848422	–	7848423	–

#### Fiber Contents:

15-20 mmHg\*: 78% Nylon, 22% Spandex.

20-30 mmHg\*: 78% Nylon, 22% Spandex.

30-40 mmHg\*: 70% Nylon, 30% Spandex.



### Knee High (Open Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518102	–	7518103	–
15-20 mmHg	<b>M</b>	7518202	–	7518203	–
15-20 mmHg	<b>L</b>	7518302	–	7518303	–
15-20 mmHg	<b>XL</b>	7518402	–	7518403	–
20-30 mmHg	<b>S</b>	7521704	–	7521705	–
20-30 mmHg	<b>M</b>	7521804	–	7521805	–
20-30 mmHg	<b>L</b>	7521904	–	7521905	–
20-30 mmHg	<b>XL</b>	7522004	–	7522005	–
20-30 mmHg	<b>LFC</b>	7848420	–	7848421	–
20-30 mmHg	<b>XLFC</b>	7848720	–	7848721	–
30-40 mmHg	<b>S</b>	7522903	–	7522904	–
30-40 mmHg	<b>M</b>	7523003	–	7523004	–
30-40 mmHg	<b>L</b>	7523103	–	7523104	–
30-40 mmHg	<b>XL</b>	7523203	–	7523204	–
30-40 mmHg	<b>LFC</b>	7848820	–	7848821	–
30-40 mmHg	<b>XLFC</b>	7848920	–	7848921	–

#### Fiber Contents:

15-20 mmHg\*: 78% Nylon, 22% Spandex.

20-30 mmHg\*: 78% Nylon, 22% Spandex.

30-40 mmHg\*: 70% Nylon, 30% Spandex.





### Petite Thigh High (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518525	–	–	–
15-20 mmHg	<b>M</b>	7518625	–	–	–
15-20 mmHg	<b>L</b>	7518725	–	–	–
15-20 mmHg	<b>XL</b>	7518825	–	–	–
20-30 mmHg	<b>S</b>	7522125	–	–	–
20-30 mmHg	<b>M</b>	7522225	–	–	–
20-30 mmHg	<b>L</b>	7522325	–	–	–
20-30 mmHg	<b>XL</b>	7522425	–	–	–
30-40 mmHg	<b>S</b>	7523325	–	–	–
30-40 mmHg	<b>M</b>	7523425	–	–	–
30-40 mmHg	<b>L</b>	7523525	–	–	–
30-40 mmHg	<b>XL</b>	7523625	–	–	–



### Petite Thigh High (Open Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518524	–	–	–
15-20 mmHg	<b>M</b>	7518624	–	–	–
15-20 mmHg	<b>L</b>	7518724	–	–	–
15-20 mmHg	<b>XL</b>	7518824	–	–	–
20-30 mmHg	<b>S</b>	7522124	–	–	–
20-30 mmHg	<b>M</b>	7522224	–	–	–
20-30 mmHg	<b>L</b>	7522324	–	–	–
20-30 mmHg	<b>XL</b>	7522424	–	–	–
30-40 mmHg	<b>S</b>	7523324	–	–	–
30-40 mmHg	<b>M</b>	7523424	–	–	–
30-40 mmHg	<b>L</b>	7523524	–	–	–
30-40 mmHg	<b>XL</b>	7523624	–	–	–



### Thigh High (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518500	7518507	7518501	7518509
15-20 mmHg	<b>M</b>	7518600	7518607	7518601	7518609
15-20 mmHg	<b>L</b>	7518700	7518707	7518701	7518709
15-20 mmHg	<b>XL</b>	7518800	7518807	7518801	7518809
20-30 mmHg	<b>S</b>	7522100	7522107	7522101	7522109
20-30 mmHg	<b>M</b>	7522200	7522207	7522201	7522209
20-30 mmHg	<b>L</b>	7522300	7522307	7522301	7522309
20-30 mmHg	<b>XL</b>	7522400	7522407	7522401	7522409
30-40 mmHg	<b>S</b>	7523300	7523307	7523301	7523309
30-40 mmHg	<b>M</b>	7523400	7523407	7523401	7523409
30-40 mmHg	<b>L</b>	7523500	7523507	7523501	7523509
30-40 mmHg	<b>XL</b>	7523600	7523607	7523601	7523609

#### Fiber Contents:

15-20 mmHg\*: 72% Nylon, 18% Spandex, 10% Silicone.

20-30 mmHg\*: 70% Nylon, 16% Spandex, 14% Silicone.

30-40 mmHg\*: 62% Nylon, 28% Spandex, 10% Silicone.

# JOBST® Opaque

## Medical Compression Stockings



### Thigh High (Open Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518502	7518508	7518503	–
15-20 mmHg	<b>M</b>	7518602	7518608	7518603	–
15-20 mmHg	<b>L</b>	7518702	7518708	7518703	–
15-20 mmHg	<b>XL</b>	7518802	7518808	7518803	–
20-30 mmHg	<b>S</b>	7522102	7522108	7522103	–
20-30 mmHg	<b>M</b>	7522202	7522208	7522203	–
20-30 mmHg	<b>L</b>	7522302	7522308	7522303	–
20-30 mmHg	<b>XL</b>	7522402	7522408	7522403	–
30-40 mmHg	<b>S</b>	7523302	7523308	7523303	–
30-40 mmHg	<b>M</b>	7523402	7523408	7523403	–
30-40 mmHg	<b>L</b>	7523502	7523508	7523503	–
30-40 mmHg	<b>XL</b>	7523602	7523608	7523603	–

#### Fiber Contents:

15-20 mmHg\*: 72% Nylon, 18% Spandex, 10% Silicone.

20-30 mmHg\*: 70% Nylon, 16% Spandex, 14% Silicone.

30-40 mmHg\*: 62% Nylon, 28% Spandex, 10% Silicone.



### Pantyhose (Closed Toe)

Compression mmHg*	Size	Natural	Honey	Black	Espresso
15-20 mmHg	<b>S</b>	7518900	–	7518901	–
15-20 mmHg	<b>M</b>	7519000	–	7519001	–
15-20 mmHg	<b>L</b>	7519100	–	7519101	–
15-20 mmHg	<b>XL</b>	7519200	–	7519201	–
20-30 mmHg	<b>S</b>	7522500	–	7522503	–
20-30 mmHg	<b>M</b>	7522600	–	7522603	–
20-30 mmHg	<b>L</b>	7522700	–	7522703	–
20-30 mmHg	<b>XL</b>	7522800	–	7522803	–
30-40 mmHg	<b>S</b>	7523700	–	7523702	–
30-40 mmHg	<b>M</b>	7523800	–	7523802	–
30-40 mmHg	<b>L</b>	7523900	–	7523902	–
30-40 mmHg	<b>XL</b>	7524000	–	7524002	–

#### Fiber Contents:

15-20 mmHg\*: 76% Nylon and 24% Spandex.

20-30 mmHg\*: 76% Nylon and 24% Spandex.

30-40 mmHg\*: 66% Nylon and 34% Spandex.

#### Care Instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size.

# JOBST<sup>®</sup> soSoft

Medical Compression Stockings



Soft, stylish and  
easy to don  
for everyday  
comfort



# JOBST® soSoft

## Medical Compression Stockings

### Features & Benefits

- Very easy to put on
- Softness and stylish pattern ensure high patient compliance
- Superior moisture management for dry and comfortable feet



### Available Colors and Styles



### Knee High (Closed Toe)

Compression mmHg*	Size	Brocade		Ribbed	
		Black	Sand	Black	Sand
15-20 mmHg	<b>S</b>	7456818	7456822	7456826	7456830
15-20 mmHg	<b>M</b>	7456819	7456823	7456827	7456831
15-20 mmHg	<b>L</b>	7456820	7456824	7456828	7456832
15-20 mmHg	<b>XL</b>	7456821	7456825	7456829	7456833
20-30 mmHg	<b>S</b>	7456834	7456838	7456842	7456846
20-30 mmHg	<b>M</b>	7456835	7456839	7456843	7456847
20-30 mmHg	<b>L</b>	7456836	7456840	7456844	7456848
20-30 mmHg	<b>XL</b>	7456837	7456841	7456845	7456849
30-40 mmHg	<b>S</b>	7456950	7457360	7457750	7458160
30-40 mmHg	<b>M</b>	7457050	7457460	7457850	7458260
30-40 mmHg	<b>L</b>	7457150	7457560	7457950	7458360
30-40 mmHg	<b>XL</b>	7457250	7457660	7458050	7458460

### Sizing Chart

Size	Ankle 	Calf 
<b>S</b>	18-21 cm	28-38 cm
<b>M</b>	21-25 cm	30-42 cm
<b>L</b>	25-29 cm	32-46 cm
<b>XL</b>	29-33 cm	34-50 cm

### Fiber Contents:

15-20 mmHg\*: 81% Nylon, 19% Spandex.  
 20-30 mmHg\*: 74% Nylon, 26% Spandex.  
 30-40 mmHg\*: 74% Nylon, 26% Spandex.

### Care Instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

Please refer to the product label and /or package insert for full instructions on the safe use of these products.

\* The mean compression for the average ankle size.

# **JOBST<sup>®</sup> forMen**

Medical Compression Stockings



**Effective leg  
therapy in a  
high-quality  
dress sock**



# JOBST<sup>®</sup> forMen

## Medical Compression Stockings

### Features & Benefits

- Fashionable ribbed design looks like a fine men's dress sock
- Reinforced heel and toe resist wear and tear
- All-day comfort band keeps socks up without binding or pinching

### Available Colors



Black



Khaki



White






Brown





Navy

### Sizing Chart

Size	Ankle 	Calf 	Thigh 
<b>S</b>	18 - 21 cm	28 - 38 cm	40 - 62 cm
<b>M</b>	21 - 25 cm	30 - 42 cm	46 - 70 cm
<b>L</b>	25 - 29 cm	32 - 46 cm	54 - 78 cm
<b>XL</b>	29 - 33 cm	34 - 50 cm	60 - 81 cm
<b>MT</b>	21 - 25 cm	30 - 42 cm	-
<b>LT</b>	25 - 29 cm	32 - 46 cm	-
<b>LFC</b>	25 - 30 cm	46 - 61 cm	-
<b>XLFC</b>	30 - 36 cm	46 - 61 cm	-

	Knee High Regular: 50 cm or less for shoe size $\geq$ 12 (CAN and US)/47 (EU)
	Knee High Tall: Greater than 50 cm for shoe size $>$ 12 (CAN and US)/47 (EU)

For Knee High, measure Ankle and Calf. For Thigh High, measure Ankle, Calf and Thigh.

### Knee High (Closed Toe)

Compression mmHg*	Size	Black	Khaki	White	Brown	Navy
15-20 mmHg	<b>S</b>	115000	115012	115008	115084	115004
15-20 mmHg	<b>M</b>	115001	115013	115009	115085	115005
15-20 mmHg	<b>L</b>	115002	115014	115010	115086	115006
15-20 mmHg	<b>XL</b>	115003	115015	115011	115087	115007
15-20 mmHg	<b>LT</b>	115068	-	-	-	-
20-30 mmHg	<b>S</b>	115088	115100	-	115096	115092
20-30 mmHg	<b>M</b>	115089	115101	-	115097	115093
20-30 mmHg	<b>L</b>	115090	115102	-	115098	115094
20-30 mmHg	<b>XL</b>	115091	115103	-	115099	115095
20-30 mmHg	<b>MT</b>	115255	-	-	-	115256
20-30 mmHg	<b>LT</b>	115260	-	-	-	115262
20-30 mmHg	<b>LFC</b>	115294	-	-	-	-
20-30 mmHg	<b>XLFC</b>	115295	-	-	-	-
30-40 mmHg	<b>S</b>	115108	115120	-	115116	115112
30-40 mmHg	<b>M</b>	115109	115121	-	115117	115113
30-40 mmHg	<b>L</b>	115110	115122	-	115118	115114
30-40 mmHg	<b>XL</b>	115111	115123	-	115119	115115
30-40 mmHg	<b>MT</b>	115257	-	-	-	115258
30-40 mmHg	<b>LT</b>	115265	-	-	-	115267
30-40 mmHg	<b>LFC</b>	115296	-	-	-	-
30-40 mmHg	<b>XLFC</b>	115297	-	-	-	-

**Fiber Contents:**

15-20 mmHg\*: 76% Nylon, 24% Spandex.  
 20-30 mmHg\*: 73% Nylon, 27% Spandex.  
 30-40 mmHg\*: 73% Nylon, 27% Spandex.

### Knee High (Open Toe)

Compression mmHg*	Size	Black	Khaki	White	Brown	Navy
20-30 mmHg	<b>S</b>	115432	-	-	-	-
20-30 mmHg	<b>M</b>	115433	-	-	-	-
20-30 mmHg	<b>L</b>	115434	-	-	-	-
20-30 mmHg	<b>XL</b>	115435	-	-	-	-
20-30 mmHg	<b>LFC</b>	115372	-	-	-	-
20-30 mmHg	<b>XLFC</b>	115373	-	-	-	-
30-40 mmHg	<b>S</b>	115452	-	-	-	-
30-40 mmHg	<b>M</b>	115453	-	-	-	-
30-40 mmHg	<b>L</b>	115454	-	-	-	-
30-40 mmHg	<b>XL</b>	115455	-	-	-	-
30-40 mmHg	<b>LFC</b>	115374	-	-	-	-
30-40 mmHg	<b>XLFC</b>	115375	-	-	-	-

**Fiber Contents:**

20-30 mmHg\*: 71% Nylon, 29% Spandex.  
 30-40 mmHg\*: 71% Nylon, 29% Spandex.

# JOBST® forMen

## Medical Compression Stockings

### Thigh High (Closed Toe)

Compression mmHg*	Size	Black	Khaki
15-20 mmHg	<b>S</b>	115516	115512
15-20 mmHg	<b>M</b>	115517	115513
15-20 mmHg	<b>L</b>	115518	115514
15-20 mmHg	<b>XL</b>	115519	115515
20-30 mmHg	<b>S</b>	115408	115400
20-30 mmHg	<b>M</b>	115409	115401
20-30 mmHg	<b>L</b>	115410	115402
20-30 mmHg	<b>XL</b>	115411	115403
30-40 mmHg	<b>S</b>	115412	115404
30-40 mmHg	<b>M</b>	115413	115405
30-40 mmHg	<b>L</b>	115414	115406
30-40 mmHg	<b>XL</b>	115415	115407

#### Fiber Contents:

15-20 mmHg\*: 73% Nylon, 20% Spandex, 7% Silicone.

20-30 mmHg\*: 71% Nylon, 25% Spandex, 4% Silicone.

30-40 mmHg\*: 71% Nylon, 25% Spandex, 4% Silicone.

**Care Instructions:** Machine washable in lukewarm water. Do not bleach. Machine dry - low heat. Do not dry clean. Do not iron. Use a delicate fabric detergent. Rinse well. Do not wring.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size.



# **JOBST<sup>®</sup> forMen Casual**

Medical Compression Stockings




**Superior comfort  
and wicking  
technology in  
a casual, versatile  
sock that fits all  
of life's occasions**



# JOBST® forMen Casual

## Medical Compression Stockings

### Features & Benefits

-  Superior moisture wicking properties to help keep feet dry and cool
- Contemporary, relaxed design that can be worn every day
- Reinforced heel and toe resist wear and tear
- All-day comfort knee band keeps socks up without binding or pinching

### Sizing chart

Size	Ankle 	Calf 
<b>S</b>	18 - 21 cm	28 - 38 cm
<b>M</b>	21 - 25 cm	30 - 42 cm
<b>L</b>	25 - 29 cm	32 - 46 cm
<b>XL</b>	29 - 33 cm	34 - 50 cm

### Available Colors



Black



Khaki

### Knee High (Closed Toe)

Compression mmHg*	Size	Black	Khaki
15-20 mmHg	<b>S</b>	7547602	7547603
15-20 mmHg	<b>M</b>	7547702	7547703
15-20 mmHg	<b>L</b>	7547802	7547803
15-20 mmHg	<b>XL</b>	7547902	7547903
20-30 mmHg	<b>S</b>	7548402	7548403
20-30 mmHg	<b>M</b>	7548502	7548503
20-30 mmHg	<b>L</b>	7548602	7548603
20-30 mmHg	<b>XL</b>	7548702	7548703
30-40 mmHg	<b>S</b>	7549202	7549203
30-40 mmHg	<b>M</b>	7549302	7549303
30-40 mmHg	<b>L</b>	7549402	7549403
30-40 mmHg	<b>XL</b>	7549502	7549503

#### Fiber Contents:

15-20 mmHg\*: 46% Nylon, 14% Spandex, 35% Polyester, 5% Wool.

20-30 mmHg\*: 47% Nylon, 16% Spandex, 33% Polyester, 4% Wool.

30-40 mmHg\*: 48% Nylon, 16% Spandex, 31% Polyester, 5% Wool

**Care Instructions:** Machine washable in lukewarm water. Do not bleach. Machine dry - low heat. Do not dry clean. Do not iron. Use a delicate fabric detergent. Rinse well. Do not wring.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

Dri-release® is an Optimer Inc. trademark.

\*The mean compression for the average ankle size.

**NEW!**

# **JOBST<sup>®</sup> forMen Ambition**

Medical Compression Stockings



**The compression stockings to which men can commit**



# JOBST® forMen Ambition

## Medical Compression Stockings

### Features & Benefits

- Special activated carbon technology integrated in foot sole to provide a unique, natural and biologically neutral odor-control-system to keep your feet fresh and cool, rechargeable every laundering
- Soft yarns and high quality cotton for easy donning and comfortable wearing
- Relief zone with flat toe seam reduces pressure points for all day wearing comfort
- New heel pocket design especially for men's anatomy for a perfect fit
- Wide top band reduces constriction and stays in place

### Available Colors



Black



Khaki



Dark Grey



Brown



Navy

### Sizing Chart

	1	2	3	4	5	6
<b>Calf</b>	28 - 36 cm	30 - 40 cm	34 - 44 cm	38 - 50 cm	44 - 56 cm	50 - 62 cm
<b>Ankle</b>	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 27 cm	27 - 30 cm	30 - 34 cm
<b>Length</b> 	Regular: 38 - 44 cm / Long: 45 - 51 cm					

### Knee High (Closed Toe)

Compression mmHg*	Size	Black		Khaki		Grey		Brown		Navy	
		Reg.	Long	Reg.	Long	Reg.	Long	Reg.	Long	Reg.	Long
15-20 mmHg	1	7766500	7766600	7766510	7766610	7766540	7766640	7766520	7766620	7766530	7766630
15-20 mmHg	2	7766501	7766601	7766511	7766611	7766541	7766641	7766521	7766621	7766531	7766631
15-20 mmHg	3	7766502	7766602	7766512	7766612	7766542	7766642	7766522	7766622	7766532	7766632
15-20 mmHg	4	7766503	7766603	7766513	7766613	7766543	7766643	7766523	7766623	7766533	7766633
15-20 mmHg	5	7766504	7766604	7766514	7766614	7766544	7766644	7766524	7766624	7766534	7766634
15-20 mmHg	6	7766505	7766605	7766515	7766615	7766545	7766645	7766525	7766625	7766535	7766635
20-30 mmHg	1	7766700	7766800	7766710	7766810	7766740	7766840	7766720	7766820	7766730	7766830
20-30 mmHg	2	7766701	7766801	7766711	7766811	7766741	7766841	7766721	7766821	7766731	7766831
20-30 mmHg	3	7766702	7766802	7766712	7766812	7766742	7766842	7766722	7766822	7766732	7766832
20-30 mmHg	4	7766703	7766803	7766713	7766813	7766743	7766843	7766723	7766823	7766733	7766833
20-30 mmHg	5	7766704	7766804	7766714	7766814	7766744	7766844	7766724	7766824	7766734	7766834
20-30 mmHg	6	7766705	7766805	7766715	7766815	7766745	7766845	7766725	7766825	7766735	7766835
30-40 mmHg	1	7766900	7767000	7766910	7767010	7766940	7767040	7766920	7767020	7766930	7767030
30-40 mmHg	2	7766901	7767001	7766911	7767011	7766941	7767041	7766921	7767021	7766931	7767031
30-40 mmHg	3	7766902	7767002	7766912	7767012	7766942	7767042	7766922	7767022	7766932	7767032
30-40 mmHg	4	7766903	7767003	7766913	7767013	7766943	7767043	7766923	7767023	7766933	7767033
30-40 mmHg	5	7766904	7767004	7766914	7767014	7766944	7767044	7766924	7767024	7766934	7767034
30-40 mmHg	6	7766905	7767005	7766915	7767015	7766945	7767045	7766925	7767025	7766935	7767035

**Fiber Contents:** Nylon, Spandex, Cotton, Polyester.

**Care Instructions:** Machine washable in lukewarm water. Do not bleach. Machine dry - low heat. Do not dry clean. Do not iron. Use a delicate fabric detergent. Rinse well. Do not wring. Do not use fabric softener.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size

# **JOBST<sup>®</sup> ActiveWear**

Medical Compression Stockings



**Effective leg  
therapy for all  
activities in  
an energizing  
athletic sock**



### Features & Benefits





- Extra cushioned foot with a reinforced heel and toe enhances wearing comfort and resists abrasion
- Seamless toe
- All-day comfort knee band keeps socks up without binding or pinching
- One unisex style to suit both men and women
- Feels like cotton and offers superior moisture management to help keep feet dry, cool and comfortable

### Knee High (Closed toe)

Compression mmHg*	Size	White	Black
15-20 mmHg	<b>S</b>	110479	110483
15-20 mmHg	<b>M</b>	110480	110484
15-20 mmHg	<b>L</b>	110481	110485
15-20 mmHg	<b>XL</b>	110482	110486
15-20 mmHg	<b>LFC</b>	110525	110531
15-20 mmHg	<b>XLFC</b>	110526	110532
20-30 mmHg	<b>S</b>	110489	110493
20-30 mmHg	<b>M</b>	110490	110494
20-30 mmHg	<b>L</b>	110491	110495
20-30 mmHg	<b>XL</b>	110492	110496
20-30 mmHg	<b>LFC</b>	110527	110533
20-30 mmHg	<b>XLFC</b>	110528	110534
30-40 mmHg	<b>S</b>	7541102	7541103
30-40 mmHg	<b>M</b>	7541202	7541203
30-40 mmHg	<b>L</b>	7541302	7541303
30-40 mmHg	<b>XL</b>	7541402	7541403
30-40 mmHg	<b>LFC</b>	110529	110535
30-40 mmHg	<b>XLFC</b>	110530	110536

### Sizing chart:

Size	Ankle 	Calf 
<b>S</b>	18-21 cm	28-38 cm
<b>M</b>	21-25 cm	30-42 cm
<b>L</b>	25-29 cm	32-46 cm
<b>XL</b>	29-33 cm	34-50 cm
<b>LFC</b>	25-30 cm	46-61 cm
<b>XLFC</b>	30-36 cm	46-61 cm

### Available Colors



White



Black



**Fiber Contents:** 15-20 mmHg\*: 60% Polyester, 20% Nylon, 11% Cotton, 9% Spandex.  
 20-30 mmHg\*: 60% Polyester, 20% Nylon, 11% Cotton, 9% Spandex.  
 30-40 mmHg\*: 58% Polyester, 20% Nylon, 10% Cotton, 12% Spandex .

**Care Instructions:** Machine washable in lukewarm water. Tumble dry in low heat. Do not bleach, dry clean, or iron.

\*The mean compression for the average ankle size

**NEW!**  
**JOBST® Sport**  
Medical Compression Stockings



**Effective leg therapy  
during all sportive  
activities**



# JOBST® Sport

## Medical Compression Stockings



### Features & Benefits

- Activated carbon yarns technology with a unique odor control system that recharges with every laundering
- Moisture control for cool and dry feet
- Padded foot helps to absorb impact during activities
- Extra-flat toe seam reduces blisters building during activities


### Knee High (Closed toe)

Compression mmHg*	Size	Light grey	Dark grey
15-20 mmHg	<b>S</b>	7528920	7528930
15-20 mmHg	<b>M</b>	7528921	7528931
15-20 mmHg	<b>L</b>	7528922	7528932
15-20 mmHg	<b>XL</b>	7528923	7528933
20-30 mmHg	<b>S</b>	7529020	7529030
20-30 mmHg	<b>M</b>	7529021	7529031
20-30 mmHg	<b>L</b>	7529022	7529032
20-30 mmHg	<b>XL</b>	7529023	7529033

### Sizing Chart

Size	Ankle 	Calf 
<b>S</b>	18-21 cm	28-38 cm
<b>M</b>	21-25 cm	30-42 cm
<b>L</b>	25-29 cm	32-46 cm
<b>XL</b>	29-33 cm	34-50 cm

### Available Colors



Light grey

Dark grey

**Fiber Contents:** 15-20 mmHg\*: 30% Polyester, 60% Nylon, 10% Spandex.  
20-30 mmHg\*: 30% Polyester, 60% Nylon, 10% Spandex.

**Care Instructions:** Machine washable in lukewarm water. Tumble dry in low heat. Do not bleach, dry clean, or iron.

\*The mean compression for the average ankle size

Please refer to the product label and/or package insert for full instructions on the safe use of these products.





# JOBST<sup>®</sup> Relief<sup>®</sup>

Medical Compression Stockings



**Compression  
therapy at an  
affordable price**



# JOBST® Relief®

## Medical Compression Stockings

### Features & Benefits

- Closed toe styles are designed with a roomy toe
- Reinforced heel and toe resist wear and tear
- Available in chap style
- Pantyhose available in open toe

### Available Colors



Beige



Black

### Sizing Chart

Size	Ankle 	Calf 	Thigh 	Hip 
<b>S</b>	18 - 21 cm	28 - 38 cm	40 - 62 cm	71 - 117 cm
<b>M</b>	21 - 25 cm	30 - 42 cm	46 - 70 cm	76 - 127 cm
<b>L</b>	25 - 29 cm	32 - 46 cm	54 - 78 cm	81 - 137 cm
<b>XL</b>	29 - 33 cm	34 - 50 cm	60 - 81 cm	102 - 165 cm
<b>LFC</b>	25 - 30 cm	46 - 61 cm	-	-
<b>XLFC</b>	30 - 36 cm	46 - 61 cm	-	-

For Knee High: Measure ankle and calf. For Thigh High: Measure ankle, calf and thigh.  
For Pantyhose and Chaps: Measure ankle, thigh and hip.

### Knee High

Compression mmHg*	Size	(Closed Toe)		(Open Toe)
		Beige	Black	Beige
15-20 mmHg	<b>S</b>	114806	114812	114800
15-20 mmHg	<b>M</b>	114807	114813	114801
15-20 mmHg	<b>L</b>	114808	114814	114802
15-20 mmHg	<b>XL</b>	114809	114815	114803
15-20 mmHg	<b>LFC</b>	114810	114816	114804
15-20 mmHg	<b>XLFC</b>	114811	114817	114805
20-30 mmHg	<b>S</b>	114620	114730	114625
20-30 mmHg	<b>M</b>	114621	114731	114626
20-30 mmHg	<b>L</b>	114622	114732	114627
20-30 mmHg	<b>XL</b>	114623	114733	114628
20-30 mmHg	<b>LFC</b>	114698	114734	114696
20-30 mmHg	<b>XLFC</b>	114624	114735	114629
30-40 mmHg	<b>S</b>	114630	114736	114635
30-40 mmHg	<b>M</b>	114631	114737	114636
30-40 mmHg	<b>L</b>	114632	114738	114637
30-40 mmHg	<b>XL</b>	114633	114739	114638
30-40 mmHg	<b>LFC</b>	114699	114740	114697
30-40 mmHg	<b>XLFC</b>	114634	114741	114639

### Fiber Contents

15-20 mmHg : CT & OT: 80% Nylon, 20% Spandex

20-30 mmHg: CT: 77% Nylon, 23% Spandex OT: 75% Nylon, 25% Spandex

30-40 mmHg: CT: 72% Nylon, 28% Spandex OT: 70% Nylon, 30% Spandex

## Thigh High With Silicone Band

Compression mmHg*	Size	(Closed Toe)		(Open Toe)
		Beige	Black	Beige
15-20 mmHg	<b>S</b>	114822	114826	114818
15-20 mmHg	<b>M</b>	114823	114827	114819
15-20 mmHg	<b>L</b>	114824	114828	114820
15-20 mmHg	<b>XL</b>	114825	114829	114821
20-30 mmHg	<b>S</b>	114208	114212	114200
20-30 mmHg	<b>M</b>	114209	114213	114201
20-30 mmHg	<b>L</b>	114210	114214	114202
20-30 mmHg	<b>XL</b>	114211	114215	114203
30-40 mmHg	<b>S</b>	114216	114220	114204
30-40 mmHg	<b>M</b>	114217	114221	114205
30-40 mmHg	<b>L</b>	114218	114222	114206
30-40 mmHg	<b>XL</b>	114219	114223	114207

### Fiber Contents

15-20 mmHg: 74% Nylon, 18% Spandex, 8% Silicone

20-30 mmHg: 70% Nylon, 23% Spandex, 7% Silicone

30-40 mmHg: 66% Nylon, 28% Spandex, 6% Silicone

## Thigh High Without Silicone Band

Compression mmHg*	Size	(Closed Toe)	(Open Toe)
		Beige	Beige
20-30 mmHg	<b>S</b>	114640	114644
20-30 mmHg	<b>M</b>	114641	114645
20-30 mmHg	<b>L</b>	114642	114646
20-30 mmHg	<b>XL</b>	114643	114647
30-40 mmHg	<b>S</b>	114648	114652
30-40 mmHg	<b>M</b>	114649	114653
30-40 mmHg	<b>L</b>	114650	114654
30-40 mmHg	<b>XL</b>	114651	114655

### Fiber Contents

20-30 mmHg:

CT: 74% Nylon, 26% Spandex

OT: 73% Nylon, 27% Spandex

30-40 mmHg:

CT: 67% Nylon, 33% Spandex

OT: 68% Nylon, 32% Spandex

# JOBST® Relief®

## Medical Compression Stockings

### Pantyhose

Compression mmHg*	Size	(Closed Toe) Beige	(Open Toe) Beige
20-30 mmHg	<b>S</b>	114656	114660
20-30 mmHg	<b>M</b>	114657	114661
20-30 mmHg	<b>L</b>	114658	114662
20-30 mmHg	<b>XL</b>	114659	114663
30-40 mmHg	<b>S</b>	114664	114668
30-40 mmHg	<b>M</b>	114665	114669
30-40 mmHg	<b>L</b>	114666	114670
30-40 mmHg	<b>XL</b>	114667	114671

#### Fiber Contents

20-30 mm Hg & 30-40 mmHg:  
CT: 75% Nylon, 25% Spandex  
OT: 78% Nylon, 22% Spandex

### Chap (Open Toe)

Compression mmHg*	Size	(Left Leg) Beige	(Right Leg) Beige	(Both Legs) Beige
20-30 mmHg	<b>S</b>	114672	114676	114680
20-30 mmHg	<b>M</b>	114673	114677	114681
20-30 mmHg	<b>L</b>	114674	114678	114682
20-30 mmHg	<b>XL</b>	114675	114679	114683
30-40 mmHg	<b>S</b>	114784	114788	114792
30-40 mmHg	<b>M</b>	114785	114789	114793
30-40 mmHg	<b>L</b>	114786	114790	114794
30-40 mmHg	<b>XL</b>	114787	114791	114795

#### Fiber Contents

20-30 mmHg & 30-40 mmHg: 78% Nylon, 22% Spandex

#### Care Instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

\*The mean compression for the average ankle size.

# JOBST® Travel Socks

Compression Stockings



For effective  
leg therapy  
that energizes  
your legs



# JOBST® Travel Socks

## Compression stockings

### Features & Benefits

- Clinically proven gradient compression helps improve blood circulation, helps prevent leg swelling, and decreases the likelihood of blood clot formation
- All-day comfort band keeps sock up
- Stylish ribbed pattern works for both business and casual wear
- Soft cotton blend provides breathable comfort
- Energizing sock design ensures comfort on arrival
- Unisex design
- Reinforced toe and heel for durability



### Ordering Information

#### Knee High (closed toe)

Compression mmHg	Size	Color	Code
15-20 mmHg	1	Black	7884407
15-20 mmHg	2	Black	7884507
15-20 mmHg	3	Black	7884607
15-20 mmHg	4	Black	7884707
15-20 mmHg	5	Black	7927207
15-20 mmHg	1	Navy	7884406
15-20 mmHg	2	Navy	7884506
15-20 mmHg	3	Navy	110355
15-20 mmHg	4	Navy	110356
15-20 mmHg	5	Navy	110357



### Sizing Chart

Size	Shoe Size		Ankle Measurement 	Calf Measurement 
	Women	Men		
1	4 1/2 - 6 1/2	3 1/2 - 5 1/2	18 - 22 cm	27 - 37 cm
2	7 - 9	6 - 8	20 - 23 cm	31 - 41 cm
3	9 1/2 - 11	8 1/2 - 10	22 - 25 cm	33 - 43 cm
4	11 1/2 - 13	10 1/2 - 12	24 - 27 cm	36 - 46 cm
5	N/A	12 1/2 - 14	25 - 29 cm	40 - 50 cm

#### Fiber Contents:

15-20 mmHg: 55% Nylon, 33% Cotton, 12% Spandex.

#### Care Instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.

# Accessories



Useful with  
compression  
garments

# Accessories

## IT STAYS®! Roll-on Body Adhesive

It Stays®! Roll-on Body Adhesive is specially designed to help keep JOBST® stockings in place. Ideal for thigh length stockings when garter belts are undesirable or silicone band is not available.



## JOLASTIC® Special Washing Solution

Washing solution for elastic garments formulated to remove soil, body acids, and skin salts quickly and easily without damage to fabric. This specially designed solution will help extend the life of all elastic garments.



## Foot Rest

Provides convenience when measuring for compression garments by elevating the patient's foot at the heel.



## Ordering Information

### It Stays®! Roll-On Body Adhesive

Code	Description	Box Contents	Case Contents
112013	2 oz. (59 mL) bottle	1 each	N/A
112014	2 oz. (59 mL) bottle	N/A	12 bottles

### Jolastic® Special Washing Solution

Code	Description	Box Contents	Case Contents
112001	12 oz. (355 mL) bottle	1 each	N/A
112003	1 qt. (946 mL) bottle	1 each	N/A
130999	4 oz. (118 mL) bottle	N/A	12 bottles

### Foot Rest

Code	Description	Box Contents	Case Contents
110912	Foot Rest	1 each	N/A



# JOBST® Donning Aids

Devices to assist patients



**Improve compliance  
and are easy to use**



# JOBST® Donning Aids

Devices to assist patients

## JOBST® Stocking Donner

This product is ideal for patients with arthritis and limited bending motion. Use with closed as well as open toe compression stockings.



## Stocking Slip

Reduces friction between the compression stocking and the skin allowing the compression stocking to slip on easily. Use with open toe compression stockings only.



## JOBST® Donning Gloves

JOBST® gloves create friction and provide grip to help donning and positioning evenly the fabric on the leg or the arm.



Code	Description	Box Contents
110913	JOBST® Stocking Donner	1 each
112018	JOBST® Stocking Slip	6 each
131202	JOBST® Donning Gloves Small	1 pair
131203	JOBST® Donning Gloves Medium	1 pair

\*Not recommended for use with Custom Seamed hosiery.





**Ask your representative or call our Customer Service  
to learn more about JOBST® Compression Stockings**

**1-877-978-5526**

**We wish you a healthy  
and comfortable day.**

