

Cover-Roll® Stretch

Gentle, secure fixation dressing



Latex-Free.
Wide-Area Fixation.

Cover-Roll® Stretch

Gentle, secure fixation dressing

Features and Benefits

When secure dressings are required, wide-area fixation now has an established place instead of surgical tapes. Wide-area fixation with Cover-Roll® Stretch combines the advantages of conventional systems because of three main reasons. First, secure and rapid fixation is possible on all, even frequently moved and highly contoured parts of the body, such as joints. Second, the entire wound dressing is protected against contamination. Third, cutting to size eliminates unnecessary waste.

- The soft, stretchable non-woven polyester material adapts well to body contours
- Provides light compression without causing constriction or congestion
- The polyacrylate adhesive holds dressing safely and reliably in place, even in presence of moisture
- Easy, secure application due to the integrated application
- Radio-transparent to eliminate the need to remove for X-ray examinations
- 100% latex-free formula to eliminate the risk of latex reactions

Indications for Use

- Retains primary dressing on small to large surgical or wound sites
- Ideal use for areas of body requiring flexibility such as elbows, knees
- Ideal base coverage for areas applying the patella/femoral taping technique



Ordering Information

Code Number	Description	Box Contents	Case Contents
45547	2" x 2 yds	1 roll	12 boxes
45548	4" x 2 yds	1 roll	12 boxes
45549	6" x 2 yds	1 roll	12 boxes
45550	8" x 2 yds	1 roll	12 boxes
45551	12" x 2 yds	1 roll	12 boxes
45552	2" x 10 yds	1 roll	12 boxes
45553	4" x 10 yds	1 roll	12 boxes
45554	6" x 10 yds	1 roll	12 boxes
45555	8" x 10 yds	1 roll	12 boxes
45556	12" x 10 yds	1 roll	12 boxes

NOTE: Please refer to the product label and/or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
 Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
 To order toll-free: BSN medical (+1) 800 552 1157